

VITA

Shade systems

Tooth shade determination

VITA shade guides



VITA Shade Guides

Beautiful teeth visible at a glance – a million times over.

The area comprising the natural tooth shades is just a small section within the entire color spectrum, however, it provides its own world of color based on the factors of lightness, chroma and hue.

The **VITA classical A1–D4®** system is a basic solution that offers 16 natural tooth shades and one that has been internationally established for decades. The **VITA SYSTEM 3D-MASTER®**, which offers 26 natural tooth shades, has also been used successfully around the world for many years and reproduces all known natural tooth shades systematically in the scientifically developed tooth color space.

VITA classical A1-D4®

VITA SYSTEM 3D-MASTER®

VITA classical A1-D4®

The basic solution with a plus.

VITA classical A1–D4® shade guide with VITA Bleached Shades

- Tooth shade determination in a single step
- 16 natural, empirically induced tooth shades
- In addition, three shades for tooth whitening using VITA Bleached Shades
- Simple shade selection thanks to quick comparison of shade samples with the natural tooth
- Possibility of sorting in accordance with visually perceived levels of lightness



VITA SYSTEM 3D-MASTER®


Linearguide, Bleachedguide and Toothguide. Efficient shade determination.



VITA Linearguide 3D-MASTER

- Tooth shade determination in two steps
- Trusted tooth shade determination thanks to classic, linear handling
- 26 natural and systematically arranged tooth shades; in addition, three shades for tooth whitening (VITA Bleached Shades)
- Clear organization of the shade samples in terms of lightness, chroma and hue
- Easy use that is self-explanatory
- Possibility of creating intermediate shades by mixing

VITA Bleachedguide 3D-MASTER®

- Instrument for determining the level of lightness for the tooth whitening process (bleaching) in just a single step
- 29 levels of lightness based on VITA SYSTEM 3D-MASTER®
- Realistic depiction of tooth whitening as well as simple communication with the patient thanks to a clear and transparent structure





VITA Toothguide 3D-MASTER® with Bleached Shades

- Tooth shade determination in three steps
- 26 natural and systematically arranged tooth shades; in addition, three shades for tooth whitening (VITA Bleached Shades)
- Clear organization of the shade samples in terms of lightness, chroma and hue
- Possibility of creating intermediate shades by mixing